

MPC MEDICAL PHYSICIST INSPECTION BEFORE FIRST USE

The following items need to be inspected by your MPC medical physicist **before** first patient use.

- *New Mammographic Unit or Processor*
- *Existing Unit or Processor Moved to a New Address*
- *Existing Unit or Processor Moved to New Location at Current Address with Major Part(s) Reassembly*
- *Automatic Exposure Control (AEC)*
 - AEC sensor replacement
 - AEC circuit board replacement
- *Collimator Replacement*
 - Entire collimator system replacement
 - Blade replacement
- *HV Generator Replacement*
- *Tube Replacement*
- *Filter Replacement*
- *Digital Software Modification*
- *Full Field Digital Mammography (FFDM) Detector - REPLACEMENT or REPAIR*
- *FFDM Monitor / Printer (NEW, ADDITIONAL, or REPLACED)*